

## Release Notes

DPNCheck® Reporter Application, version 2.3

**Note: For a more detailed description overview of new features, please download the DPNCheck Reporter User Manual for version 2.3 located at [www.dpncheck/Reporter](http://www.dpncheck/Reporter).**

1. Cloud Reporting: Introduction of a new feature, which, when enabled, allows for sending of de-identified Reporter study data to the DPNCheck Cloud for viewing in aggregate with a Dashboard. A new “Cloud” tab is added to settings to support the feature, and a SQLite database is saved locally (stores de-identified data only). An active internet connection is required for use of this feature. This feature is disabled by default.
2. EHR Integration: Introduction of a new feature, which when enabled, allows for sending an HL7 export message directly to an HL7 MLLP receiver over TCP/IP to a defined IP address and port number. The User Interface is updated to indicate status of HL7 transmission. A new “EHR” tab is added to settings to support the feature. This feature is disabled by default.
3. Normal Limits: Addition of a new Normal Limits option that only requires entry of Birth Date for limit calculation. Use of normal limits remains optional and under user control.
4. Multiple enhancements for usability, including:
  - NPI number displayed on Report tab.
  - Inclusion/exclusion of Practice Information entry on report tab controlled by new setting.
  - Modified display of Patient to group required fields at the top; height field only displayed if required for Normal Limit calculation.
  - When selected from Physician List, fields associated with Physician may not be edited.
  - Improvements in messages to users.
5. Load-time Improvements: Ability to control time period for studies loaded on start-up in Settings.
6. Logging: Logging to support new functionality, integration with Microsoft Eventlog, default logging contains no PII/PHI.
7. Removal of select feature support:
  - XML export format
  - Weight entry
  - Use of First and/or Last Name as part of study filename; neither may be used
  - Spanish language
  - Auto-ID population
  - Normal Limits confirmation no longer required during installation process
  - Datasheet
  - Use of Export from within Study Search
  - Day-Month date format selection options
  - Merge into a Completed Study (Now Merge is only allowed for New Studies, only applies if combining multiple tests with same Patient ID into the same study)
  - Local storage of Monograph and User Manual pdfs; accessible through website
8. Branding updates: Logo, color




# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3

### Reporter 2.1.1.4

Help menu label change from “Documentation” to “Resources”, goes to same website destination.  
Modified About screen, new branding. Reporter no longer has a version of the User Manual distributed with the app; all resources are on website.




# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3




### Reporter 2.3

New tabs for EHR and Cloud, Application tab removed.



### Reporter 2.1.1.4


Application settings removed (Language, Auto population ID, File name).



# Summary of Primary Screen Changes


## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3





### Reporter 2.1.1.4

*Date Range option removed.*



*Enable Practice Entry Setting added.*




# Summary of Primary Screen Changes


## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3



No changes



### Reporter 2.1.1.4



*EMR Export Format XML removed, Auto Export relocated to EHR tab.*




# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3

New tab for EHR settings, Automatic Export relocated to this tab, new setting to support EHR Export to IP Address destination.

### Reporter 2.1.1.4



When Export to EHR enabled, Study List shows status of EHR Export.





# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3


### Reporter 2.3

### Reporter 2.1.1.4

*No changes.*



*New tab for Cloud settings, multiple settings to support Cloud feature.*



# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3

### Reporter 2.1.1.4

*New age-only normal limit option.*

Normal Limits Definition

Normal limits for DPNCheck are described in the DPNCheck Normative Data monograph, which can be found at [www.dpncheck.com/resources](http://www.dpncheck.com/resources).

- Use age and height adjusted normal limits (from monograph).
- Use age-only adjusted normal limits (from monograph).
- Use standardized normal limits (5  $\mu$ V, 41 m/s).
- Use custom normal limits.
- Do not use normal limits.

Continue

Normal Limits Definition

Optional normal limits for DPNCheck are described in the following monograph: [DPNCheck Normative Data, Collection, Analysis and Recommended Normal Limits](#)

- Use age and height adjusted normal limits (from monograph).
- Use standardized normal limits (5  $\mu$ V, 41 m/s).
- Use custom normal limits.
- Do not use normal limits.


Continue

# Summary of Primary Screen Changes


DPNCheck® Reporter Application, version 2.3

## Reporter 2.3

*Patient Tab: Reorganization of fields with required fields at the top and designated with \*, Height field only displayed if required for normal limits calculation.*



## Reporter 2.1.1.4




# Summary of Primary Screen Changes

## DPNCheck® Reporter Application, version 2.3

### Reporter 2.3

*Report Tab: NPI field displayed. Practice field display (Practice, Address, City/State/Zip) may be disabled in Settings. When selected from Physician List, fields associated with Physician may not be edited.*



### Reporter 2.1.1.4

